

# EQUIPMENT RENTAL REQUEST: LYMPHASYS-9



Email, fax or post this completed form to Medi-Rent and a member of our team will call you to complete your order. Keep a copy for yourself for referrers recommended settings. Alternatively, order online at [www.medirent.com.au](http://www.medirent.com.au)

## 1. ORDER DETAILS

Only reqd for insurance

Required Delivery Date

-----

Clinician

-----

Signature

-----

Provider No.

-----

Phone No.

-----

## 2. PATIENT DETAILS

Patient Name

-----

Address

-----

-----

Home phone

-----

Mobile

-----

## 3. LYMPHASYS-9

### Treatment Plan

Your clinician must complete this section for your prescribed settings. (Refer overleaf for setting options)

Limb length

-----

Limb circumference at widest point

-----

### PRE-TREATMENT MODE A

Pressure:

-----

Session Duration:

-----

### TREATMENT MODE B

Pressure:

-----

Speed:

-----

Session Duration:

-----

Times a day:

-----

## 4. BILLING INFORMATION

Who is paying for the equipment?

Hospital

Client

Insurance (please fill form below)

**Pricing: \$120 per month**, includes garments  
 Please call us to discuss terms, conditions, and delivery / collection options.

## 5. INSURANCE INFORMATION

If you are claiming this equipment through your private health insurance, please make a second copy of this form to send to your health fund.

Insurance Company

-----

Claim No.

-----

Phone No.

-----

Contact Name

-----

Address

-----

-----

Email

-----

Request Approval Received?

Yes

No

Received By:  Email

Fax

Letter

**PH: 1300 726 666**

# LYMPHASYS-9

## WHAT IS SIPC?

Sequential Intermittent Pneumatic Compression (SIPC) therapy involves the use of external pressure to promote the circulation of blood and lymphatic fluid to and from the extremities of the body. Each chamber inflates and deflates sequentially starting from the extremity and moving towards the centre of the body.


It is beneficial in the prevention and treatment of many circulatory problems such as lymphoedema, chronic venous insufficiency & venous ulcers, as well as sports recovery.

Medi-Rent offers a Rent-Try-Buy scheme where we will refund up to 3 months rental off the total purchase price.




OPTIONS	
Pre-Treatment Mode - A (lymphoedema only)	The proximal area is cleared to prepare for distal oedema to drain efficiently
Normal Circulation Mode - B	Fluid is pushed sequentially starting distally towards to torso
Pressure Range	10 - 180mmHg
Session Duration	5 - 90 minutes
Speed Range	1 - 6 mL/s
Mode	A-D

### A MODE PRE-TREATMENT



### B MODE TREATMENT



*IPC reduces venous stasis and increases flow velocity in the deep veins resulting in favourable haemodynamic changes such as decreased venous pressure and interstitial oedema.*

- Comerota (2011)

