

RENTAL REQUEST FORM: LX9 (Sequential intermittent pneumatic compression)

Email, fax or post this completed form to Medi-Rent and a member of our team will call you to complete your order. Keep a copy for yourself for referrers recommended settings. Alternatively, order online at www.medirent.com.au

1. ORDER DETAILS

Required Delivery Date

Clinician

Clinician's address

Phone No.

2. PATIENT DETAILS

Patient Name

Address

Home phone

Mobile

3. BILLING INFORMATION

Who is paying for the equipment?

Hospital

Local Health District

Client

LX9 TREATMENT PLAN -

Your clinician must complete this section for your prescribed settings. (Refer overleaf for setting and sizing options)

MEASUREMENTS -

Only required for lower limbs as the arm and abdominal garments are a universal size.

Garment/s required: Leg Arm Abdomen

Leg garment size if known: -----

Otherwise complete the relevant sections

Leg length - from groin to ankle:

Left ----- Right -----

Leg circumference - top of thigh:

Left ----- Right -----

Leg circumference - ankle:

Left ----- Right -----

PRE-TREATMENT MODE A -

Lymphoedema only. See overleaf for more detail

Pressure: -----

Session duration: -----

TREATMENT MODE B

Pressure: -----

Speed: -----

Session duration: -----

Times a day: -----

Pricing: \$130 per month, includes required leg or arm garments. Plus delivery. Additional garments will incur an extra fee.

Please call us to discuss terms, conditions, and delivery / collection options.


PH: 1300 726 666

LX9 - SETTINGS

OPTIONS


Pre-Treatment Mode - A (Lymphoedema only)	The proximal area is cleared to prepare for distal oedema to drain efficiently
Normal Circulation Mode - B	Fluid is pushed sequentially starting distally towards to torso
Pressure Range	10 - 180mmHg
Session Duration	5 - 90 minutes
Speed Range	1 - 6 mL/s

LX9 ABDOMINAL GARMENT



One size


LX9 ARM GARMENT



	A	B
One size	50cm	88cm
with extension zipper	60cm	88cm

*error range : ±2cm

LX9 LEG GARMENT



	Thigh (+extension)	Ankle (+extension)	Length
M	50cm (60cm)	36cm (46cm)	67cm
L	62cm (72cm)	37cm (44cm)	74cm
XL	74cm (84cm)	43cm (50cm)	79cm
XXL	79cm (89cm)	46cm (53cm)	92cm

PH: 1300 726 666