

# RENTAL REQUEST FORM: LX9 (Sequential intermittent pneumatic compression)

Email, fax or post this completed form to Medi-Rent and call us with payment details to complete your order. Alternatively, order online at [www.medirent.com.au](http://www.medirent.com.au)

## 1. ORDER DETAILS

Date Requested

-----

Required Delivery Date

-----

Clinician

-----

1. Courier to:  (Fees apply)

Clinician's address

-----

Phone No.

-----

## 2. PATIENT DETAILS

Patient

-----

DOB

-----

2. Courier to:  (Fees apply)

Patient's Address

-----

Email

-----

Home phone

-----

Mobile

-----

or 3. Home set-up Sydney Metro (Premium service fees apply)

## 3. BILLING INFORMATION

Who is paying?

Hospital -----

Client

## LX9 REQUIREMENTS

Clinicians must complete this section for prescribed settings and garments

### TREATMENT OPTIONS

Mode - A Pre-Treatment	Pre-clearing of proximal area in prep for efficient fluid drainage	Pressure Range 10 - 140mmHg
Mode - B Classic Sequential	Fluid is pushed sequentially distally to proximally	
Mode - C Distal oedema	Four pumps on distal chamber before inflating sequentially	Session Duration 5 - 90 minutes, or continuous
Mode - D Wave	Garment fills two chambers at any one time in a wave formation, followed by classic sequential	Cycle Speed 1 - 6 (slowest - fastest)

### PRE-TREATMENT MODE A

Pressure: -----

Session duration: -----

Speed: -----

TREATMENT MODE B  MODE C  MODE D

Pressure: -----

Session duration: -----

Speed: -----

### GARMENTS *To determine sizing see overleaf*

Garment/s reqd:  Legs  Arm  Abdomen

Leg garment size: -----

Number of extensions: ----- pair/s

## MEASUREMENTS

Measure your patient in relevant areas to determine sizing

1. Leg length - from groin to ankle:

Left ..... Right .....

2. Thigh circumference (widest part):

Left ..... Right .....

3. Calf circumference:


Left ..... Right .....

4. Ankle circumference:

Left ..... Right .....

5. Trunk circumference (widest part): .....

6. Affected arm circumference (widest part): .....



## LEG GARMENT

	Length	Ankle	Calf	Thigh
M	68cm	30cm	37cm	47cm
+1 ext	68cm	37cm	45cm	57cm
+2 ext	68cm	44cm	53cm	67cm
+3 ext	68cm	51cm	61cm	77cm

L	74cm	36cm	45cm	59cm
+1 ext	74cm	43cm	53cm	69cm
+2 ext	74cm	50cm	61cm	79cm
+3 ext	74cm	57cm	69cm	89cm

XL	78cm	41cm	50cm	70cm
+1 ext	78cm	48cm	58cm	80cm
+2 ext	78cm	55cm	66cm	90cm
+3 ext	78cm	62cm	74cm	100cm


XXL	93cm	42cm	56cm	73cm
+1 ext	93cm	49cm	64cm	83cm
+2 ext	93cm	56cm	72cm	93cm
+3 ext	93cm	63cm	80cm	103cm

## TRUNK GARMENT

One size with 3 built-in zippers

	On widest fitting	+ 1 pair extensions	+ 2 pairs extensions	+ 3 pairs extensions
B	133cm	153cm	173cm	193cm
C	74cm	84cm	94cm	104cm

## ARM GARMENT



	A = Upper limb circumference
One size fits all garment	50cm
+ 1 extension	60cm
+ 2 extensions	70cm