



ABOUT US

Medi-Rent offers a range of specialist medical equipment for rent to medical institutions, practitioners and your patients. Our 25 years of experience has given us a comprehensive understanding of our medical equipment, so we know we can genuinely deliver the very best information and service to our clients.

We believe everyone has a right to the best rehabilitation.

Recovery from a hospital visit, surgery or a course of therapy usually doesn't end when the patient returns home. But buying expensive equipment to fully recover at home is not a viable option for many. That is why our business is focussed on providing equipment solutions that are traditionally known as 'expensive equipment' at competitive rental rates. We offer direct monthly billing and a rent-try-buy scheme for some products.

Looking after our customers is key.

We pride ourselves on our detailed knowledge of our products including their operation, application and service requirements. Ordering with Medi-Rent is quick and easy with knowledgeable customer service representatives ready to take your call. We deliver Australia-wide and can provide set-up and instruction.

Commitment to clinical efficacy.

Our equipment is hand-picked to fit with our company values and code of ethics. It's important to us that our products are grounded by good science so that we can stand by their claims to provide improved health outcomes. We continue to stay abreast of the latest clinical research in the industry to ensure we are providing credible solutions. Medi-Rent is committed to providing education and in-servicing to healthcare providers on the use and application of our systems.

SPECIFICATIONS

Variable pressure	10-140mmHg
Modes of Compression	Pre-treatment massage & therapy modes
Power	AC 220-240V, 50/60Hz
Time Range	5 - 90 minutes, or continuous
Product Weight	2kg
Product Size	260 x 160 x 120mm

Easy delivery and set-up options at affordable prices.



CONTACT US

Medi-Rent Pty Ltd

Ph: 1300 726 666
Intl: + 612 9700 1744
F: + 612 9316 9998
E: info@medirent.com.au

www.medirent.com.au
www.facebook.com/pages/Medirent



VASCULAR COMPRESSION PUMPS

Sequential Intermittent Pneumatic Compression





VASCULAR COMPRESSION PUMPS

Intermittent Pneumatic Compression mimics the body's natural lymphatic drainage.

- ✓ Heals skin ulceration
- ✓ Softens secondary skin changes
- ✓ Reduces oedema volume

- ✓ LIGHT
- ✓ PORTABLE
- ✓ COMFORTABLE TO WEAR
- ✓ EASY TO PUT ON



LX9 MODES

- A mode - Pre-therapy, lymph node clearance
- B mode - classic sequential compression
- C mode - distal intensive setting
- D mode - wave and sequential sequence

SIPC

Sequential Intermittent Pneumatic Compression Therapy (SIPC) is the process of applying external pressure to the limbs to promote the circulation of blood and lymphatic fluid to and from the extremities of the body. It may be beneficial for the management of the following conditions:

- Lymphoedema
- Chronic venous insufficiency
- Venous ulceration
- Lipoedema
- Diabetic foot ulceration

HOW IT WORKS

Intermittent pneumatic compression mimics the body's manual lymphatic drainage. By applying pressure to the surface of the limb it promotes the circulation of blood/fluid throughout the body. Each chamber inflates and deflates sequentially starting from the extremity and moving towards the centre of the body.

BENEFITS

- Can be more cost effective long term than intensive lymphoedema massage or bandaging.
- Can be used in the comfort of the home.
- A good alternative to compression stockings or bandaging when those products aren't suitable.
- Garments are included in the price of the rental with the device, unlike competitors who charge additional fees.

Our compression pump has four chambers which is most cost effective and portable for home use without compromising on clinical outcomes.

RENT TRY BUY SCHEME
If you decide to purchase, we refund up to **3 MONTHS** of paid rental off the purchase price.

LX9 FEATURES

- Variable pressure 10-140mmHg
- Continuous mode - 80mmHg max pressure
- Therapy compliance timer (on request)
- Lockable programs (on request)
- Easy fit garments - leg, arm or trunk
- Various sizes and extension zippers
- Four treatment modes
- Six cycle speeds
- Sanitary, washable boot-liners (optional extra)



If you'd like to review the evidence for SIPC, please contact us for a copy of our Clinical Guidelines.

